

Lactotransferrin

Catalog No: #AP79883



Package Size: #AP79883-1 50ug #AP79883-2 100ug #AP79883-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lactotransferrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Lactoferrin
Accession No.	Uniprot ID:P24627
Uniprot	P24627
GeneID	280846;
SDS-PAGE MW	22.44kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:40-243aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LTF

Note: This product is for in vitro research use only