

Protein arginine N-methyltransferase 1

Catalog No: #AP79884



Package Size: #AP79884-1 50ug #AP79884-2 100ug #AP79884-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein arginine N-methyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Histone-arginine N-methyltransferase PRMT1, Interferon receptor 1-bound protein 4, HMT2, HRMT1L2, IR1B4
Accession No.	Uniprot ID:Q99873
Uniprot	Q99873
GeneID	3276;
SDS-PAGE MW	18.7kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:41-210aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)gene_full_name:PRMT1

Note: This product is for in vitro research use only