

Chromodomain-helicase-DNA-binding protein 8

Catalog No: #AP79890



Package Size: #AP79890-1 50ug #AP79890-2 100ug #AP79890-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chromodomain-helicase-DNA-binding protein 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CHD-8,ATP-dependent helicase CHD8,Helicase with SNF2 domain 1,HELSNF1,KIAA1564
Accession No.	Uniprot ID:Q9HCK8
Uniprot	Q9HCK8
GeneID	57680;
SDS-PAGE MW	29.7kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:451-720aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CHD8

Note: This product is for in vitro research use only