

Macrosialin

Catalog No: #AP79891



Package Size: #AP79891-1 50ug #AP79891-2 100ug #AP79891-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Macrosialin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Gp110,CD68
Accession No.	Uniprot ID:P34810
Uniprot	P34810
GeneID	968;
SDS-PAGE MW	12.87kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:48-164aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD68

Note: This product is for in vitro research use only