

# Non-histone chromosomal protein HMG-17

Catalog No: #AP79893



Package Size: #AP79893-1 50ug #AP79893-2 100ug #AP79893-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Non-histone chromosomal protein HMG-17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	High mobility group nucleosome-binding domain-containing protein 2,HMG17
Accession No.	Uniprot ID:P80272
Uniprot	P80272
GeneID	100524500;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HMGN2

Note: This product is for in vitro research use only