

# Insulin-like growth factor I

Catalog No: #AP79894

Package Size: #AP79894-1 50ug #AP79894-2 100ug #AP79894-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Insulin-like growth factor I
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IGF-I,Somatomedin
Accession No.	Uniprot ID:P07455
Uniprot	P07455
GeneID	281239;
SDS-PAGE MW	11.55kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:50-154aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IGF1

Note: This product is for in vitro research use only