

LIM domain and actin-binding protein 1

Catalog No: #AP79895



Package Size: #AP79895-1 50ug #AP79895-2 100ug #AP79895-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LIM domain and actin-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Epithelial protein lost in neoplasm,EPLIN,SREBP3,PP624
Accession No.	Uniprot ID:Q9UHB6
Uniprot	Q9UHB6
GeneID	51474;
SDS-PAGE MW	14.41kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:51-181aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LIMA1

Note: This product is for in vitro research use only