

NADPH oxidase 3

Catalog No: #AP79897



Package Size: #AP79897-1 50ug #AP79897-2 100ug #AP79897-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NADPH oxidase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Mitogenic oxidase 2,MOX-2,gp91phox homolog 3,GP91-3,MOX2
Accession No.	Uniprot ID:Q9HBY0
Uniprot	Q9HBY0
GeneID	50508;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NOX3

Note: This product is for in vitro research use only