

Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit alpha isoform



Catalog No: #AP79898

Orders: order@signalwayantibody.com

Package Size: #AP79898-1 50ug #AP79898-2 100ug #AP79898-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit alpha isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PI3-kinase subunit alpha,Phosphatidylinositol 4,5-bisphosphate 3-kinase 110 kDa catalytic subunit alpha,PtdIns-3-kinase subunit p110-alpha,Phosphoinositide-3-kinase catalytic alpha polypeptide,Serine/threonine protein kinase PIK3CA
Accession No.	Uniprot ID:P42336
Uniprot	P42336
GeneID	5290;
SDS-PAGE MW	23.21kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:520-730aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PIK3CA

Note: This product is for in vitro research use only