

Laminin subunit beta-1

Catalog No: #AP79899



Package Size: #AP79899-1 50ug #AP79899-2 100ug #AP79899-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Laminin subunit beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Laminin B1 chain,Laminin-1 subunit beta,Laminin-10 subunit beta,Laminin-12 subunit beta,Laminin-2 subunit beta,Laminin-6 subunit beta,Laminin-8 subunit beta,Lamb-1,Lamb1-1
Accession No.	Uniprot ID:P02469
Uniprot	P02469
GeneID	16777;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lamb1

Note: This product is for in vitro research use only