

Fibrinogen gamma chain

Catalog No: #AP79903



Package Size: #AP79903-1 50ug #AP79903-2 100ug #AP79903-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibrinogen gamma chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PRO2061
Accession No.	Uniprot ID:P02679
Uniprot	P02679
GeneID	2266;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGG

Note: This product is for in vitro research use only