

# Myogenic factor 6

Catalog No: #AP79904



Package Size: #AP79904-1 50ug #AP79904-2 100ug #AP79904-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Myogenic factor 6                                      |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Pig  |
| Other Names        | Myf-6  |
| Accession No.      | Uniprot ID:Q3YFL6                                      |
| Uniprot            | Q3YFL6   |
| GeneID             | 397005;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MYF6

Note: This product is for in vitro research use only