

Protein 4.1

Catalog No: #AP79905



Package Size: #AP79905-1 50ug #AP79905-2 100ug #AP79905-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein 4.1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	P4.1,4.1R,Band 4.1,EPB4.1,E41P
Accession No.	Uniprot ID:P11171
Uniprot	P11171
GeneID	2035;
SDS-PAGE MW	27.94kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:611-864aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EPB41

Note: This product is for in vitro research use only