

# Protein 4.1

Catalog No: #AP79906



Package Size: #AP79906-1 50ug #AP79906-2 100ug #AP79906-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein 4.1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	P4.1,4.1R,Band 4.1,Epb4.1,Kiaa4056
Accession No.	Uniprot ID:P48193
Uniprot	P48193
GeneID	269587;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Epb41

Note: This product is for in vitro research use only