

Antizyme inhibitor 2

Catalog No: #AP79909



Package Size: #AP79909-1 50ug #AP79909-2 100ug #AP79909-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Antizyme inhibitor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Arginine decarboxylase,ADC,ARGDC,ODC-paralogue,ODC-p,Ornithine decarboxylase-like protein,KIAA1945,ODCP,AZIN2
Accession No.	Uniprot ID:Q96A70
Uniprot	Q96A70
GeneID	113451;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADC

Note: This product is for in vitro research use only