

Tumor necrosis factor receptor superfamily member 17

Catalog No: #AP79910

Package Size: #AP79910-1 50ug #AP79910-2 100ug #AP79910-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B-cell maturation protein,CD269,Bcm,Bcma
Accession No.	Uniprot ID:O88472
Uniprot	O88472
GeneID	21935;
SDS-PAGE MW	12.65kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:71-185aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnfrsf17

Note: This product is for in vitro research use only