

# Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 1 protein

Catalog No: #AP79911

Package Size: #AP79911-1 50ug #AP79911-2 100ug #AP79911-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 1 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	P-Rex1,Kiaa1415
Accession No.	Uniprot ID:Q69ZK0
Uniprot	Q69ZK0
GeneID	277360;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Prex1

Note: This product is for in vitro research use only