

Delta-sarcoglycan

Catalog No: #AP79913



Package Size: #AP79913-1 50ug #AP79913-2 100ug #AP79913-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Delta-sarcoglycan
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Delta-SG,35 kDa dystrophin-associated glycoprotein,35DAG
Accession No.	Uniprot ID:Q92629
Uniprot	Q92629
GeneID	6444;
SDS-PAGE MW	20.46kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:72-257aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SGCD

Note: This product is for in vitro research use only