

Serine/threonine-protein kinase Sgk2

Catalog No: #AP79914



Package Size: #AP79914-1 50ug #AP79914-2 100ug #AP79914-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase Sgk2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Serum/glucocorticoid-regulated kinase 2
Accession No.	Uniprot ID:Q8R4U9
Uniprot	Q8R4U9
GeneID	171497;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sgk2

Note: This product is for in vitro research use only