

Ewing's tumor-associated antigen 1 homolog

Catalog No: #AP79915



Package Size: #AP79915-1 50ug #AP79915-2 100ug #AP79915-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ewing's tumor-associated antigen 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q5SVT3
Uniprot	Q5SVT3
GeneID	68145;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Etaa1

Note: This product is for in vitro research use only