

# Tumor necrosis factor receptor superfamily member 17

Catalog No: #AP79916

Package Size: #AP79916-1 50ug #AP79916-2 100ug #AP79916-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor receptor superfamily member 17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	B-cell maturation protein,CD269,BCM,BCMA
Accession No.	Uniprot ID:Q02223
Uniprot	Q02223
GeneID	608;
SDS-PAGE MW	11.77kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:78-184aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TNFRSF17

Note: This product is for in vitro research use only