

## Poly [ADP-ribose] polymerase 1

Catalog No: #AP79917



Package Size: #AP79917-1 50ug #AP79917-2 100ug #AP79917-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Poly [ADP-ribose] polymerase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PARP-1,ADP-ribosyltransferase diphtheria toxin-like 1,ARTD1,NAD(+) ADP-ribosyltransferase 1,ADPRT 1,Poly[ADP-ribose] synthase 1,ADPRT,PPOL
Accession No.	Uniprot ID:P09874
Uniprot	P09874
GeneID	142;
SDS-PAGE MW	16.5kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:81-230aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PARP1

Note: This product is for in vitro research use only