

Glucosylceramidase

Catalog No: #AP79918



Package Size: #AP79918-1 50ug #AP79918-2 100ug #AP79918-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glucosylceramidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Acid beta-glucosidase,Alglucerase,Beta-glucocerebrosidase,Beta-GC,D-glucosyl-N-acylsphingosine glucohydrolase,Imiglucerase,GC,GLUC
Accession No.	Uniprot ID:P04062
Uniprot	P04062
GeneID	2629;
SDS-PAGE MW	35.2kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:81-400aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GBA

Note: This product is for in vitro research use only