

EMILIN-1

Catalog No: #AP79919



Package Size: #AP79919-1 50ug #AP79919-2 100ug #AP79919-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EMILIN-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Elastin microfibril interface-located protein 1,Elastin microfibril interfacier 1,EMI
Accession No.	Uniprot ID:Q9Y6C2
Uniprot	Q9Y6C2
GeneID	11117;
SDS-PAGE MW	23.54kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:87-300aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EMILIN1

Note: This product is for in vitro research use only