

V-set domain-containing T-cell activation inhibitor 1

Catalog No: #AP79924



Package Size: #AP79924-1 50ug #AP79924-2 100ug #AP79924-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	V-set domain-containing T-cell activation inhibitor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	B7 homolog 4,B7-H4,B7h.5,Immune costimulatory protein B7-H4,Protein B7S1,T-cell costimulatory molecule B7x,B7H4,UNQ659/PRO1291
Accession No.	Uniprot ID:Q7Z7D3
Uniprot	Q7Z7D3
GeneID	79679;
SDS-PAGE MW	17.16kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:90-245aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:VTCN1

Note: This product is for in vitro research use only