

CAAX prenyl protease 1 homolog

Catalog No: #AP79925



Package Size: #AP79925-1 50ug #AP79925-2 100ug #AP79925-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CAAX prenyl protease 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Farnesylated proteins-converting enzyme 1,FACE-1,Prenyl protein-specific endoprotease 1,Zinc metalloproteinase Ste24 homolog,Face1
Accession No.	Uniprot ID:Q80W54
Uniprot	Q80W54
GeneID	230709;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Zmpste24

Note: This product is for in vitro research use only