

Multidrug resistance protein 1

Catalog No: #AP79926



Package Size: #AP79926-1 50ug #AP79926-2 100ug #AP79926-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Multidrug resistance protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ATP-binding cassette sub-family B member 1,P-glycoprotein 1,Abcb1b,Mdr1,Mdr1b,Pgy1,Pgy1-1
Accession No.	Uniprot ID:P06795
Uniprot	P06795
GeneID	18669;
SDS-PAGE MW	31.13kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:994-1276aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abcb1

Note: This product is for in vitro research use only