

## GFP Monoclonal Antibody(Mix)

Catalog No: #44046

Package Size: #44046-1 50ul #44046-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	GFP Monoclonal Antibody(Mix)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Mix
Isotype	IgG1
Applications	WB IP
Specificity	The antibody detectsThe antibody detects RFP, YFP, CFP.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	GFP
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:500-1:2000

Immunoprecipitation: 1:200


## Images



Western blot analysis of GFP transfected HeLa, diluted at

1) 1:5000

2) 1:10000



- 1) Input (control)
- 2) IP products, antibody dilution 1:200

Note: This product is for in vitro research use only