

HERV-K_7p22.1 provirus ancestral Gag polyprotein

Catalog No: #AP79931



Package Size: #AP79931-1 50ug #AP79931-2 100ug #AP79931-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HERV-K_7p22.1 provirus ancestral Gag polyprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EnvK2 protein,Envelope polyprotein,HERV-K(C7) envelope protein,HERV-K(HML-2.HOM) envelope protein,HERV-K108 envelope protein
Accession No.	Uniprot ID:Q69384
Uniprot	Q69384
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:CMNPSEMQRKAPPPRRRRHRNRStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM

Imidazole,10%glycerol(PH8.0)gene_full_name:ERVK6

Note: This product is for in vitro research use only