

# Solute carrier family 2, facilitated glucose transporter member 4

Catalog No: #AP79936

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP79936-1 50ug #AP79936-2 100ug #AP79936-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Solute carrier family 2, facilitated glucose transporter member 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Glucose transporter type 4, insulin-responsive, GLUT-4, GLUT4, B146
Accession No.	Uniprot ID:Q9XST2
Uniprot	Q9XST2
GeneID	403673;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: SLC2A4

Note: This product is for in vitro research use only