

Prokineticin-1

Catalog No: #AP79937



Package Size: #AP79937-1 50ug #AP79937-2 100ug #AP79937-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prokineticin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Endocrine-gland-derived vascular endothelial growth factor,EG-VEGF,PK1
Accession No.	Uniprot ID:Q14A28
Uniprot	Q14A28
GeneID	246691;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prok1

Note: This product is for in vitro research use only