

GFP Monoclonal Antibody(Plant Specific)

Catalog No: #44050

Package Size: #44050-1 50ul #44050-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	GFP Monoclonal Antibody(Plant Specific)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Plant Specific
Isotype	IgG1
Applications	WB IP
Specificity	The antibody detects GFP, EGFP and GFP, EGFP tag fusion proteins in Plants.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	GFP
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Immunoprecipitation: 1:200

Images



- 1) Input: Hela Cell Lysate
- 2) IP product: IP dilute 1:200

Note: This product is for in vitro research use only