

N-terminal propeptide of Collagen alpha-1(I) chain

Catalog No: #AP79945



Package Size: #AP79945-1 50ug #AP79945-2 100ug #AP79945-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | N-terminal propeptide of Collagen alpha-1(I) chain |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PINP

Note: This product is for in vitro research use only