

7-methylguanosine phosphate-specific 5'-nucleotidase

Catalog No: #AP79954



Package Size: #AP79954-1 50ug #AP79954-2 100ug #AP79954-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	7-methylguanosine phosphate-specific 5'-nucleotidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	7-methylguanosine nucleotidase, Cytosolic 5'-nucleotidase 3B, Cytosolic 5'-nucleotidase III-like protein, cN-III-like protein, N(7)-methylguanylate 5'-phosphatase, Nt5c3l
Accession No.	Uniprot ID:Q3UFY7
Uniprot	Q3UFY7
GeneID	68106;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nt5c3b

Note: This product is for in vitro research use only