

Adenosine receptor A1

Catalog No: #AP79957

Package Size: #AP79957-1 50ug #AP79957-2 100ug #AP79957-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenosine receptor A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	RDC7
Accession No.	Uniprot ID:P11616
Uniprot	P11616
GeneID	403961;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ADORA1

Note: This product is for in vitro research use only