

N-terminal propeptide of Collagen alpha-1(II) chain

Catalog No: #AP79962



Package Size: #AP79962-1 50ug #AP79962-2 100ug #AP79962-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	N-terminal propeptide of Collagen alpha-1(II) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PIINP

Note: This product is for in vitro research use only