

14-3-3 protein zeta/delta

Catalog No: #AP79965



Package Size: #AP79965-1 50ug #AP79965-2 100ug #AP79965-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	14-3-3 protein zeta/delta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Factor activating exoenzyme S,FAS,Protein kinase C inhibitor protein 1,KCIP-1
Accession No.	Uniprot ID:P63103
Uniprot	P63103
GeneID	287022;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:YWHAZ

Note: This product is for in vitro research use only