

Alpha-1-antiproteinase

Catalog No: #AP79972



Package Size: #AP79972-1 50ug #AP79972-2 100ug #AP79972-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-1-antiproteinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alpha-1-antitrypsin,Alpha-1-proteinase inhibitor,Serpin A1
Accession No.	Uniprot ID:P17475
Uniprot	P17475
GeneID	24648;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Serpina1

Note: This product is for in vitro research use only