

## Myc-Tag Monoclonal Antibody(Mix)

Catalog No: #44054

Package Size: #44054-1 50ul #44054-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Myc-Tag Monoclonal Antibody(Mix)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Mix
Isotype	IgG1
Applications	WB
Specificity	The antibody can detect C-Myc fusion proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Myc-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:500-1:2000

## Images



Western blot analysis of His-MBP-C-Myc recombinant protein,  
diluted at  
1) 1:5000  
2) 1:10000

Note: This product is for in vitro research use only