

# 3-hydroxyanthranilate 3,4-dioxygenase

Catalog No: #AP79980



Package Size: #AP79980-1 50ug #AP79980-2 100ug #AP79980-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	3-hydroxyanthranilate 3,4-dioxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	3-hydroxyanthranilate oxygenase,3-HAO,3-hydroxyanthranilic acid dioxygenase,HAD
Accession No.	Uniprot ID:P46953
Uniprot	P46953
GeneID	56823;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Haao

Note: This product is for in vitro research use only