

HLA class II histocompatibility antigen, DRB1-9 beta chain

Catalog No: #AP79987

Package Size: #AP79987-1 50ug #AP79987-2 100ug #AP79987-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | HLA class II histocompatibility antigen, DRB1-9 beta chain |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | MHC class II antigen DRB1*9,DR-9 |
| Accession No. | Uniprot ID:Q9TQE0 |
| Uniprot | Q9TQE0 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HLA-DRB1

Note: This product is for in vitro research use only