

## V5-Tag Monoclonal Antibody(Mix)

Catalog No: #44055

Package Size: #44055-1 50ul #44055-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	V5-Tag Monoclonal Antibody(Mix)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Mix
Isotype	IgG1
Applications	WB
Specificity	The antibody detects C-terminal, internal, and N-terminal V5-tag fusion proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	V5-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:500-1:2000

## Images



Western blot analysis of His-MBP-V5 recombinant protein, diluted at 1:10000.

Note: This product is for in vitro research use only