

5-hydroxytryptamine receptor 2A

Catalog No: #AP79997



Package Size: #AP79997-1 50ug #AP79997-2 100ug #AP79997-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	5-hydroxytryptamine receptor 2A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	5-HT-2,Serotonin receptor 2A,Htr2
Accession No.	Uniprot ID:P14842
Uniprot	P14842
GeneID	29595;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Htr2a

Note: This product is for in vitro research use only