

HA-Tag Monoclonal Antibody(2C1)

Catalog No: #44057

Package Size: #44057-1 50ul #44057-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HA-Tag Monoclonal Antibody(2C1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2C1
Isotype	IgG1
Applications	WB
Specificity	The antibody detects C-terminal, internal, and N-terminal HA-tag fusion proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	HA-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of Recombinant protein, diluted at

1) 1:5000

2) 1:10000

Note: This product is for in vitro research use only