

Protein 4.1

Catalog No: #AP80012



Package Size: #AP80012-1 50ug #AP80012-2 100ug #AP80012-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein 4.1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	P4.1,4.1R,Band 4.1
Accession No.	Uniprot ID:Q6Q7P4
Uniprot	Q6Q7P4
GeneID	442955;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EPB41

Note: This product is for in vitro research use only