

5'-nucleotidase

Catalog No: #AP80016



Package Size: #AP80016-1 50ug #AP80016-2 100ug #AP80016-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	5'-nucleotidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	5'-NT,Ecto-5'-nucleotidase,CD73,Nt5,Nte
Accession No.	Uniprot ID:P21588
Uniprot	P21588
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nt5e

Note: This product is for in vitro research use only