

5-hydroxytryptamine receptor 1B

Catalog No: #AP80020



Package Size: #AP80020-1 50ug #AP80020-2 100ug #AP80020-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	5-hydroxytryptamine receptor 1B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	5-HT-1B,5-HT1B,5-HTR1B,5-HT1D subtype beta,Serotonin receptor 1B
Accession No.	Uniprot ID:P79250
Uniprot	P79250
GeneID	403741;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HTR1B

Note: This product is for in vitro research use only