

5-hydroxytryptamine receptor 2C

Catalog No: #AP80040



Package Size: #AP80040-1 50ug #AP80040-2 100ug #AP80040-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	5-hydroxytryptamine receptor 2C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	5-HT-2C,5-hydroxytryptamine receptor 1C,5-HT-1C,Serotonin receptor 2C,HTR1C
Accession No.	Uniprot ID:P28335
Uniprot	P28335
GeneID	3358;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HTR2C

Note: This product is for in vitro research use only