

Cytosolic 5'-nucleotidase III-like protein

Catalog No: #AP80042



Package Size: #AP80042-1 50ug #AP80042-2 100ug #AP80042-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cytosolic 5'-nucleotidase III-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	cN-III-like protein
Accession No.	Uniprot ID:Q3UFY7
Uniprot	Q3UFY7
GeneID	68106;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nt5c3l

Note: This product is for in vitro research use only