

14-3-3 protein zeta/delta

Catalog No: #AP80052



Package Size: #AP80052-1 50ug #AP80052-2 100ug #AP80052-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	14-3-3 protein zeta/delta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Mitochondrial import stimulation factor S1 subunit,Protein kinase C inhibitor protein 1,KCIP-1,Msfs1
Accession No.	Uniprot ID:P63102
Uniprot	P63102
GeneID	25578;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ywhaz

Note: This product is for in vitro research use only